




SURGIX

DANBARROW™ MICROSURGICAL SET

Catalog of Instruments



Dan Barrow, MD

"The foundational set for microsurgical procedures"

Table of Contents

1	SET LAYOUT
2-3	SCISSORS
4-6	FORCEPS
7-8	DISSECTORS AND HOOKS
9	SUCTIONS
10	INSTRUMENT TRAYS
11	SYSTEM SUMMARY

DanBarrow™ Basic Microsurgery Set

Designed for the microsurgeon with precision and efficiency in mind, the DanBarrow™ Basic Microsurgery Set brings together the most commonly used instruments for complex microsurgery into a streamlined, all-in-one system. By consolidating the instruments that are typically spread across multiple trays, we've optimized the workflow for both the surgeon and the operating room staff. This innovation reduces unnecessary tray handling, enhances turnover speed, and boosts overall OR efficiency. Each instrument is master-crafted, and the tray is highlighted with the single-shaft dissectors that are finely balanced to the surgeon's natural hand movement and optimal center of gravity via our TrueBalance™ technology. The Dan Barrow Microscissors are designed with the right mix of tension that is firm enough for superior dissection yet carries a controlled softness for precise delicate cutting.




TrueBalance™ Technology

Developed in collaboration with Ray White and exclusively licensed to Surgix, DanBarrow Instruments' TrueBalance™ (Patent Pending) technology is engineered to elevate surgical performance. By precisely defining the ergonomic balance point of each individual instrument and counterbalancing it at the handle, TrueBalance™ delivers exceptional tactile precision and optimal weight distribution, reducing surgeon fatigue while enhancing control and accuracy.



The Foundation for Microsurgical Success

- Intuitive design
- Enhanced tactile feel
- Integrated tray design
- Custom microscissors



G2DB7000 - DBX - DanBarrow™ Basic Microsurgery Set



G2DB7011 - STRAIGHT MICROSCISSORS, BAYONET ROUND HANDLE

7" Overall Length, 3" Working Length, 10mm Tip



G2DB7012 - CURVED UP MICROSCISSORS, BAYONET ROUND HANDLE

7" Overall Length, 3" Working Length, 10mm Tip



G2DB7015 - STRAIGHT MICROSCISSORS, BAYONET ROUND HANDLE

8" Overall Length, 4" Working Length, 12mm Tip



G2DB7016 CURVED MICROSCISSORS, BAYONET ROUND HANDLE

8" Overall Length, 4" Working Length, 12mm Tip




G2DB7019 - STRAIGHT MICROSCISSORS, ROUND HANDLE

7" Overall Length, 3" Working Length, 12mm Tip



G2DB7020 - DBX STRAIGHT MICROSCISSORS, ROUND HANDLE, CURVED

7" Overall Length, 3" Working Length, 12mm Tip



G2DB7009 - DBX BARROW VASCULAR SCISSORS CERAMIC, LARGE

7" Overall Length, 5.5" Working Length



G2DB7010 - DBX BARROW VASCULAR SCISSORS CERAMIC, MEDIUM

5" Overall Length, 3.75" Working Length



G2DB7023 - DBX SUTURE TYING FCPS, STR, 0.6mm TITANIUM

6" Overall Length, 1.25" Working Length



G2DB7024- DBX SUTURE TYING FCPS, CVD, 0.6mm TITANIUM

6" Overall Length, 1.25" Working Length




G2DB7025 - DBX SUTURE TYING FCPS, STR, 0.3mm TITANIUM

6" Overall Length, 1.25" Working Length



G2DB7026 - DBX SUTURE TYING FCPS, CVD, 0.3mm TITANIUM

6" Overall Length, 1.25" Working Length




G2DB7027 - DBX BARROW FCP, 6", TITANIUM

6" Overall Length, 1.25" Working Length




G2DB7028 - DBX FCP, 1MM RING TIP, DIAMOND DUST, TITANIUM

6" Overall Length, 1.25" Working Length



G2DB7030 - DBX MICRO FCP, ROUND BAYO, STR, 0.5mm TIP, MEDIUM, CERAMIC
7.5" Overall Length, 3" Working Length



G2DB7031 - DBX MICRO FCP, ROUND BAYO, STR, 0.5MM TIP, LONG, CERAMIC
8.25" Overall Length, 3.75" Working Length



G2DB7032- DBX MICRO FCP, ROUND BAYO, STR, 1MM TIP, CERAMIC
7.5" Overall Length, 3" Working Length



G2DB7033 - DBX MICRO FCP, ROUND BAYO, STR, 1MM TIP, CERAMIC
8.25" Overall Length, 3.75" Working Length



G2DB7029- DBX MICRO FCP, 1X2 TEETH, CERAMIC
7.5" Overall Length, 2.5" Working Length



G2DB7006 - DBX RING HANDLE MICRO CUP FCP, STR, 2MM CUP, CERAMIC
7" Overall Length, 5" Working Length




G2DB7037 - DBX BARROW TUMOR CUP FCP, STR, 2MM CUP, w/TOOTH, CERAMIC
10" Overall Length, 6.5" Working Length



G2DB7001 - DBX MICRO SPATULA, SMALL, TRUEBALANCE CERAMIC

7.5" Overall Length, 3.25" Working Length



G2DB7002 - DBX MICRO SPATULA, MEDIUM, TRUEBALANCE CERAMIC

7.5" Overall Length, 3.25" Working Length



G2DB7003 - DBX MICRO SPATULA, LARGE, TRUEBALANCE CERAMIC

7.5" Overall Length, 3.25" Working Length



G2DB7004- DBX MICRO BALL DISSECTOR, 4mm, 45d ANG TRUEBALANCE CERAMIC

7.5" Overall Length, 3.25" Working Length



G2DB7005 - DBX MICRO BALL DISSECTOR, 8mm, 45d ANG TRUEBALANCE CERAMIC

7.5" Overall Length, 3.25" Working Length



G2DB7035 - DBX NERVE HOOK, TRUEBALANCE CERAMIC

7.5" Overall Length, 3.25" Working Length



G2DB7036 - DBX DOUBLE ENDED TUMOR CURETTE, CERAMIC

6.5" Overall Length, 1.75" Working Length




G2DB7034 - DBX BALL NERVE HOOK, TRUEBALANCE CERAMIC

7.5" Overall Length, 3.25" Working Length




G2DB7008 - DBX SUTURE PUSHER, TRUEBALANCE CERAMIC

7.5" Overall Length, 3.25" Working Length



G2DB7045 - DBX TUMOR HOOK, STANDARD, TRUEBALANCE CERAMIC

8" Overall Length, 4" Working Length



G2DB7046 - DBX TUMOR HOOK, EXTRA FINE, TRUEBALANCE CERAMIC

8" Overall Length, 4" Working Length



G2DB7040- DBX ATRAUMATIC CONTROLLED SUCTION, 4FR, CERAMIC
7" Overall Length, 4" Working Length





G2DB7041- DBX ATRAUMATIC CONTROLLED SUCTION, 5FR, CERAMIC
7" Overall Length, 4" Working Length



G2DB7042- DBX ATRAUMATIC CONTROLLED SUCTION, 6FR, CERAMIC
7" Overall Length, 4" Working Length



DanBarrow Suction's unique design combines the benefit of an atraumatic tip with the superior control of the teardrop control port.



Thoughtfully designed
right down to the trays.

INSTRUMENT TRAYS

DanBarrow™ Basic Microsurgery

Set Summary

SURGIX Item#	DESCRIPTION	Overall Length	Working Length
G2DB7000	DBX DANBARROW™ BASIC MICROSURGERY SET		
G2DB7001	DBX MICRO SPATULA, SMALL, TRUE BALANCE CERAMIC	7.5"	3.25"
G2DB7002	DBX MICRO SPATULA, MEDIUM, TRUE BALANCE CERAMIC	7.5"	3.25"
G2DB7003	DBX MICRO SPATULA, LARGE, TRUE BALANCE CERAMIC	7.5"	3.25"
G2DB7004	DBX MICRO BALL DISSECTOR, 4MM, 45d ANG TRUE BALANCE CERAMIC	7.5"	3.25"
G2DB7005	DBX MICRO BALL DISSECTOR, 8MM, 45d ANG TRUE BALANCE CERAMIC	7.5"	3.25"
G2DB7006	DBX RING HANDLE MICRO CUP FCP, STR, 2MM CUP, SHORT, CERAMIC	7"	5"
G2DB7008	DBX SUTURE PUSHER, TRUE BALANCE CERAMIC	7.5"	3.25"
G2DB7009	DBX BARROW VASCULAR SCISSORS CERAMIC, LARGE	7"	5.5"
G2DB7010	DBX BARROW VASCULAR SCISSORS CERAMIC, MEDIUM	5"	3.75"
G2DB7011	DBX MICROSCISSOR ROUND BAYO, REG BLADE, STR, SHORT, CERAMIC	7"	3"
G2DB7012	DBX MICROSCISSOR ROUND BAYO, REG BLADE, CVD, SHORT, CERAMIC	7"	3"
G2DB7015	DBX MICROSCISSOR ROUND BAYO, REG BLADE, STR, LONG, CERAMIC	8"	4"
G2DB7016	DBX MICROSCISSOR ROUND BAYO, REG BLADE, CVD, LONG, CERAMIC	8"	4"
G2DB7019	DBX MICROSCISSORS, ROUND, REG BLADE, STR, MEDIUM, CERAMIC	7"	3"
G2DB7020	DBX MICROSCISSORS, ROUND, REG BLADE, CVD, MEDIUM, CERAMIC	7"	3"
G2DB7023	DBX SUTURE TYING FCPS, STR, 0.6MM TITANIUM	6"	1.25"
G2DB7024	DBX SUTURE TYING FCPS, CVD, 0.6MM TITANIUM	6"	1.25"
G2DB7025	DBX SUTURE TYING FCPS, STR, 0.3MM TITANIUM	6"	1.25"
G2DB7026	DBX SUTURE TYING FCPS, CVD, 0.3MM TITANIUM	6"	1.25"
G2DB7027	DBX BARROW FCP, 6", TITANIUM	6"	1.25"
G2DB7028	DBX FCP, 1MM RING, DIAMOND DUST, TITANIUM	6"	1.25"
G2DB7029	DBX MICRO FCP, 1X2 TEETH, LONG, CERAMIC	7.5"	2.5"
G2DB7030	DBX MICRO FCP, ROUND BAYO, STR, 0.5MM TIP, MEDIUM, CERAMIC	7.5"	3"
G2DB7031	DBX MICRO FCP, ROUND BAYO, STR, 0.5MM TIP, LONG, CERAMIC	8.25"	3.75"
G2DB7032	DBX MICRO FCP, ROUND BAYO, STR, 1MM TIP, MEDIUM, CERAMIC	7.5"	3"
G2DB7033	DBX MICRO FCP, ROUND BAYO, STR, 1MM TIP, LONG, CERAMIC	8.25"	3.75"

DanBarrow™ Basic Microsurgery Set

Summary Continued-

G2DB7034	DBX BALL NERVE HOOK, TRUEBALANCE™ CERAMIC	7.5"	3.25"
G2DB7035	DBX NERVE HOOK, TRUEBALANCE™ CERAMIC	7.5"	3.25"
G2DB7036	DBX DOUBLE ENDED TUMOR CURETTE, CERAMIC	6.5"	1.75"
G2DB7037	DBX BARROW TUMOR CUP FCP, STR, 2MM CUP, LONG, w/TOOTH, CERAMIC	10"	6.5"
G2DB7040	DBX ATRAUMATIC CONTROLLED SUCTION, 4FR, CERAMIC	7"	4"
G2DB7041	DBX ATRAUMATIC CONTROLLED SUCTION, 5FR, CERAMIC	7"	4"
G2DB7042	DBX ATRAUMATIC CONTROLLED SUCTION, 6FR, CERAMIC	7"	4"
G2DB7045	DBX TUMOR HOOK - STANDARD, TRUEBALANCE CERAMIC	8"	4"
G2DB7046	DBX TUMOR HOOK - EXTRA FINE, TRUEBALANCE CERAMIC	8"	4"
G2DB7050	DBX STERILIZATION CASE WITH BUILT IN DISSECTION RACK		



DanBarrow™ Basic Microsurgery Set

The goal was to design a foundational microsurgery set that serves as an extension of the surgeon's skill, experience, and training—providing both basic and advanced instruments to maximize the efficiency of the microscope and the precision of microsurgical techniques.

*TrueBalance™ technology was developed in conjunction with Ray White and licensed exclusively to Surgix Medical, LLC. *Patent Pending.



Email. info@surgixmedical.com
Phone. (888) 247-0477
Fax. (346) 510-4823
Web. SurgixMedical.com

Request a demo.

