



FUSIONMAX® ELITE TABLE



ONE TABLE

SPINE | ORTHO | PAIN | GENERAL SURGERY



UNOBSTRUCTED C-ARM & O-ARM® ACCESS

OPEN & FLAT OPTIONS

Constructed from carbon fiber, the Open and Flat tabletops provide full radiolucency



Open design allows unrestricted imaging access

OPEN TOP

SPINE PROCEDURES



QUICK-CHANGE™ TOP

PATENT PENDING



In under a minute, a single nurse can change the tabletop on the **FusionMAX® Elite Table** from Open to Flat expanding table use beyond spine procedures for greater versatility.

FLAT TOP

ORTHO, PAIN & GENERAL SURGERY PROCEDURES



FUSIONMAX® ELITE TABLE SPECS

FUSIONMAX® BASE

Table Top to Floor Height	22-38 in. (56-96 cm)
Trendelenburg/Reverse Trend	12°
Lateral Tilt	20°
Base Width	21 in. (53 cm)
Base Length	80.5 in. (204 cm)



FUSIONMAX® QUICK-CHANGE™ FLAT TOP

Table Top Length	80.5 in. (204 cm)
Table Top Width	19.75 in. (50 cm)
Patient Weight Capacity	600 lbs. (272 kg)



FUSIONMAX® SPINAL SURGERY OPEN TOP

Table Top Length	80.5 in. (204 cm)
Table Top Width	19.25 in. (49 cm)
Patient Weight Capacity	600 lbs. (272 kg)



Email. info@surgixmedical.com

Phone. (888) 247-0477

Fax. (346) 510-4823

Web. SurgixMedical.com

Request a demo.

*O-ARM® is a registered trademark of Medtronic, Inc. Ziehm Vision RFD is a registered trademark of Ziehm Imaging.

